




# COVID-19 Antigen Rapid Test Kit

**15 MINUTES**  
DETECTION TIME

 European CE certified



## Benefits Of antigen rapid test

- **Fast**
- **No Queues, No Wait time** for results
- **Low Cost**
- **Early Detection & Elimination**
- Effort less, and **accessible** to everyone nationwide, instantly
- Viable **precautionary** step
- Objective, accurate results
- Read without instrument





# Detection Period 15 MINUTES


Excellent Performance compared with PCR

For Professional Use

Works With Nasal/Throat Swab


Tests for Antigen Simultaneously

POCT Screening



## Diagnostic Process

COVID-19 Antigen Rapid Test Kit uses advanced lateral flow technology in a sandwich design for qualitative detection of nucleocapsid protein from SARS-CoV-2. COVID-19 Antigen Rapid Test Kit, provides objective results in 15 minutes, allowing for testing of patients suspected of COVID-19 in near-patient testing environments.







## Packing Information

<b>Product Name</b>	Antigen Rapid Test Kit
<b>Product Code</b>	601460
<b>Format</b>	25T/kit
<b>Dim.</b>	17 × 12.5 × 7cm
<b>N.W.</b>	0.3kg
<b>Box Dim.</b>	70.3 × 51.3 × 36.3cm
<b>Box N.W.</b>	15.8kg
<b>Number of Box</b>	50kits

## Technical Comparison

	Fluorescence PCR	ELISA	Colloidal Gold method	Colloidal Gold method
<b>Detection substance</b>	Nucleic acid	Antibody	Antibody	Antigen
<b>Type of sample</b>	nasopharyngeal swabs, sputum, alveolar lavage fluid	Serum/plasma	serum/plasma/whole blood	Nasal/Throat Swab
<b>Time to get result</b>	2 hrs	1 hrs	Fastest 10 minutes	15 minutes
<b>Instrument needed or not</b>	Yes	Yes	No	No
<b>Laboratory requirement</b>	Confirming diagnosis	Nucleic acid negative sample reviewing or high-volume sample detection	Nucleic acid negative sample reviewing or basic hospital sample testing	Early Diagnosis

## Ordering Info

Code	Product Name	Format	Methodology
601460	Covid-19 Antigen Rapid Test Kit	25T/Kit	Colloidal Gold

### Beijing Kewei Clinical Diagnostic Reagent Inc.

**Address :** No.7, Yan Qi He Xi Yi Rd., Huai Rou District, Beijing, 101407, China

**Tel :** 86 10 68863176 / 68863259

**Website :** <http://www.keweidiagnostic.com>

**Email :** [kewei@keweidiagnostic.com](mailto:kewei@keweidiagnostic.com)

